



# High Performance Cooled CCD Camera System ASCENT<sup>®</sup> A2000

The Ascent A2000 has a 2-megapixel interline transfer sensor with high quantum efficiency and twice the dynamic range of similar Sony sensors. The camera is available with a monochrome sensor (A2000) or color sensor with a Bayer mosaic (A2000C). Low noise and small pixels are ideal for OEM applications, biological sciences, and astronomy.

- Chemiluminescence
- Fluorescence
- Astronomy
- Colorimetry



Imaging Area of CCD

- 1600 x 1200 array, 7.4 x 7.4 micron pixels
- 32 Mbyte SDRAM image buffer
- Programmable 16-bit digitization speeds up to 16 Mpixels/sec
- Single or dual channel readout
- Video focus mode
- USB 2.0 interface: no plug in cards or external controllers
- Programmable, intelligent cooling to 35°C below ambient
- Binning up to 8 Horizontal x 1200 Vertical
- Subarray readout and fast sequencing modes
- Programmable offset and gain
- External triggering and strobe controls
- ActiveX drivers included with every system
- Field upgradeable firmware
- Fused silica window
- Optional C-mount, Nikon F-mount, or 2" slip fit adapter
- Optional 6-position 1"/25mm or 8-position 1.25"/31mm filter wheel
- Single 6V supply
- Compact enclosure: 20 oz. (0.57 kg)
- Programmable status indicators



## CCD SPECIFICATIONS (from CCD manufacturer)

CCD	Kodak KAI-2020
Array Size (pixels)	1600 x 1200
Pixel Size	7.4 x 7.4 microns
Imaging Area	11.8 x 8.9 mm (105 mm <sup>2</sup> )
Imaging Diagonal	14.8 mm
Video Imager Size	0.93"
Linear Full Well (typical)	40K electrons
Dynamic Range	74 dB
QE at 400 nm (mono)	47%
Peak QE (480 nm)	56%
Anti-blooming (nominal)	300X



Part Numbers	
A2000 (mono)	A1D-02020MS-FS
A2000C (color)	A1D-02020CS-FS

For complete CCD specifications, including cosmetic grading, see data sheet from manufacturer.



151 N. Sunrise Ave. Ste 902  
Roseville CA 95661 USA  
tel 916 218 7450  
fax 916 218 7451  
www.ccd.com

PC Interface	USB 2.0
USB2 Cable	Std.: 5m. Extensions: 5 meters between hubs; 5 hubs maximum (max. total of 30m) Wide variety of extenders available, including fiber optics to 10 km.
Digital Resolution	16 bits: 3.3 and 8 Mpixels/sec (single channel); 16 Mpixels/sec (dual 8 MHz) (software selectable)
System Noise (typical)	6 e <sup>-</sup> RMS at 3.3 MHz; 7 e <sup>-</sup> RMS at 8/16 MHz
Pixel Binning	1 x 1 to 8 x 1200 on-chip
Exposure Time	Minimum 100 microseconds; max.183 minutes
Image Sequencing	1 to 65535 image sequences under software control
Frame Sizes	Full frame, subframe, focus mode
Cooling (typical)	Thermoelectric cooler. Maximum forced air cooling 35°C below ambient temperature.
Dark Current (typical)	0.01 eps (-10°C)
Temperature Stability	± 0.1°C
Camera Head Size	Aluminum. 3.2" x 4.7" x 1.3" (8.1 x 11.9 x 3.3 cm) Weight: 1.25 lb. (0.57 kg)
Mounting	1.5" x 2.5" bolt pattern, 6-32 thread. Optional C-mount (1" 32 tpi thread), Nikon F-mount, or 2" slip-fit adapters.
Back Focal Distance	Standard: 0.32" (0.81 cm). [optical]
Operating Environment	-30°C to 35°C. Relative humidity: 10 to 90% non-condensing.
Op.Sys.Support	Windows, Linux, Mac OSX
Power	20W maximum power with internal shutter open and cooling maximum. AC/DC "brick" supply with int'l AC input plug (100-240V, 50-60 Hz). Alternate 6V input from user's source.
Remote Triggering	LVTTL input allows exposure to start within 25 microseconds of rising edge of trigger

### CCD SENSITIVITY

